CURTISS - WRIGHT

R-net DIAGNOSTIC GUIDE

If the Programmer displays **Unrecognized Error**, this is probably because the Programmer software needs upgrading to include the Help Text for that Error Code. Please refer to the Programmer's documentation for details of how to upgrade to the latest software.

[DUTCH](#)

[FRENCH](#)

[GERMAN](#)

[ITALIAN](#)

[NORWEGIAN](#)

[SPANISH](#)

[SWEDISH](#)

[DANISH](#)

[SLOVAK](#)

[PORTUGUESE](#)

[FINNISH](#)

[CZECH](#)

[GREEK](#)

Diagnostic Warning

Diagnostics should only be conducted by electronic service professionals with in-depth knowledge of PG Drives Technology electronic controllers. An incorrect or badly effected repair could result in an unsafe set-up of a vehicle.

PG Drives Technology accepts no liability for losses of any kind arising from an incorrect or badly effected repair.

PG Drives Technology accept no liability for losses of any kind arising from unauthorised opening, adjustments or modifications to any component of a control system.

Module Error

The module reporting the error may require repair. Check all cables and connectors to the module, and cycle the power. If the error persists, contact your service center.

Memory Error

One of the modules in the system may require repair. Check all cables and connectors to the modules, and cycle the power.

If the error persists, disconnect each module in turn, each time cycling the power, to see if the problem can be isolated to a specific module.

If the error is corrected by disconnecting a particular module, then contact your service center about that module.

If the error cannot be corrected by disconnecting a particular module, use the R-net PC Programmer to re-program the system. This should be carried out with either: the current, specific program for the wheelchair; or the manufacturer's original program.

If, after re-programming, the error persists, contact your service center.

Warning: Programming should only be conducted by healthcare professionals with in-depth knowledge of PGDT control systems. Incorrect programming could result in an unsafe set-up of the wheelchair for a user. PGDT accept no responsibility for losses of any kind if the programming of control system is altered from the factory pre-set values.

Power Module Memory Error

An error has been detected within the Power Module. Check all cable and connectors to the Power Module, and cycle the power.

If the error persists, use the R-net PC Programmer to re-program the system. This should be carried out with either: the current, specific program for the wheelchair; or the manufacturer's original program.

If, after re-programming, the error persists, contact your service center.

Warning: Programming should only be conducted by healthcare professionals with in-depth knowledge of PGDT control systems. Incorrect programming could result in an unsafe set-up of the wheelchair for a user. PGDT accept no responsibility for losses of any kind if the programming of control system is altered from the factory pre-set values.

Cal Error

A calibration error has been detected. It may be possible to clear this error by performing the Joystick Calibration procedure in On-Board Programming.

If the error persists, contact your service center.

Joystick Error

The module reporting the error has detected a fault in the joystick that is connected to it. In the case of the Joystick Module, it may be permissible, depending on your service arrangement with the wheelchair manufacturer, to replace the joystick within the Joystick Module. If this is permitted, ensure you follow the correct anti-static and calibration procedures.

If you are not permitted to replace the joystick, then contact your service center.

In the case of specialty control interface modules, such as the Omni or PGDT I/O Module, reporting the error, then the Specialty Input Device should be replaced.

If the error persists, contact your service center.

Gone to Sleep

The system has gone to sleep because of a period of inactivity that is greater than the programmed value of Sleep Timer. To awake the system, cycle the power.

M1 Brake Error

The Power Module has detected that the brake connected via its M1 output has become disconnected. Whether the M1 output is for the left or the right brake depends on the configuration of the wheelchair and the programming of the Power Module.

Ensure there are no freewheel levers disengaged, also check all cables and connections between the Power Module M1 connector and the brake, then cycle the power.

If the error persists, contact your service center.

M2 Brake Error

The Power Module has detected that the brake connected via its M2 output has become disconnected. Whether the M1 output is for the left or the right brake depends on the configuration of the wheelchair and the programming of the Power Module.

Ensure there are no freewheel levers disengaged, also check all cables and connections between the Power Module M2 connector and the brake, then cycle the power.

If the error persists, contact your service center.

High Battery

The module reporting the error has detected that the battery voltage has risen above 35V. The most common reasons for this are: overcharging of the battery; or bad connections between the battery and the Power Module.

Ensure the charger is operating correctly and check all cables and connections between the battery and the Power Module.

If the error persists, contact your service center.

Charging

The module reporting this condition has detected that a charger is connected to it. This will be via the off-board charging socket on the Joystick Module or will be because an on-board charger has been connected to the AC line.

To clear the system, disconnect an off-board charger from the wheelchair or the on-board charger from the AC, and cycle the power.

If the condition persists, contact your service center.

Low Battery

The module reporting the error has detected a low battery voltage. The level of this voltage will differ between modules, as each can tolerate different minimum levels.

Charge the batteries. If the problem persists, check the condition of the batteries and all cables and connectors to them, and re-charge.

If the error persists, contact your service center.

Bad Settings

This occurs when the system detects incorrect or invalid programming data.

Check the programmed settings of the system, particularly those in the Profiles and Configuration groups.

If you cannot identify any settings that may be incorrect, save all the current settings in a file, then re-program the system using the manufacturer's default settings and cycle the power. If the error persists, the Power Module may be defective and you should contact your service center.

If the above procedure cleared the error, then use the saved file and program each group of parameters individually. When the error re-occurs, it will be a parameter in the last programmed group that caused it. You should therefore closely scrutinise all the settings in that group.

M1 Motor Error

The Power Module has detected that the motor connected via its M1 output has become disconnected. Whether the M1 output is for the left or the right motor depends on the configuration of the wheelchair and the programming of the Power Module.

Check all cables and connections between the Power Module M1 connector and the motor, then cycle the power.

If the error persists, contact your service center.

M2 Motor Error

The Power Module has detected that the motor connected via its M2 output has become disconnected. Whether the M2 output is for the left or the right motor depends on the configuration of the wheelchair and the programming of the Power Module.

Check all cables and connections between the Power Module M2 connector and the motor, then cycle the power.

If the error persists, contact your service center.

Overtemp. (Lamps)

The Intelligent Seating/lighting Module (ISM) has detected that its internal lighting circuit(s) have become too hot. The ISM will cease powering the lamp that is drawing too much current.

Check the wiring and connectors to the lamp in question, allow the ISM to cool down and cycle the power.

If the error persists, contact your service center.

Overtemp. (Acts)

The Intelligent Seating/lighting Module (ISM) has detected that its internal actuator circuit(s) have become too hot. The ISM will cease powering the actuator channel that is drawing too much current.

Check the wiring and connectors to the channel in question, allow the ISM to cool down and cycle the power.

If the error persists, check the actuator mechanism is in good working order. If there are any obstructions or there is too much friction in the system, then more heat will be generated within the ISM.

If the error persists, contact your service center.

Brake Lamp Short

The Intelligent Seating/lighting Module (ISM) has detected an over-current condition in its brake lamp output(s).

Check the wiring and connectors to the brake lamp(s) and cycle the power.

If the error persists, contact your service center.

Left Lamp Short

The Intelligent Seating/lighting Module (ISM) has detected an over-current condition on its left-hand driving lamp output(s).

Check the wiring and connectors to the left-hand driving lamps(s) and cycle the power.

If the error persists, contact your service center.

L Ind Lamp Short

The Intelligent Seating/lighting Module (ISM) has detected an over-current condition on its left-hand indicator output.

Check the wiring and connectors to the left-hand indicator(s) and cycle the power.

If the error persists, contact your service center.

Ind Lamp Failed

The Intelligent Seating/lighting Module (ISM) has detected an error in the wheelchair's indicator circuits.

If you have access to the four digit error codes, then 882C suggests that the connection to both indicators on a particular side has been lost, and 883C suggests that just a single bulb has failed.

Check each indicator bulb, and the cables and connectors to each of the indicators, and cycle the power.

If the error persists, contact your service center.

Right Lamp Short

The Intelligent Seating/lighting Module (ISM) has detected an over-current condition on its right-hand driving lamp output(s).

Check the wiring and connectors to the right-hand driving lamps(s) and cycle the power.

If the error persists, contact your service center.

R Ind Lamp Short

The Intelligent Seating/lighting Module (ISM) has detected an over-current condition in its right-hand indicator output.

Check the wiring and connectors to the right-hand indicator(s).

If the error persists, contact your service center.

Over-current

The Intelligent Seating/lighting Module (ISM) has detected an over-current condition in the wheelchairs actuator circuit(s).

Check the wiring and connectors to all the wheelchair's actuators and cycle the power.

If the problem persists, check the actuator mechanisms are in good working order. If there are any obstructions or there is too much friction in the system, then an excessive amount of current may be drawn. Check also that any end-stop switches are operating correctly and terminate the power to an actuator.

If the problem persists, contact your service center.

Bad Cable

The system has detected a bad connection in one of the cables that connect the modules. Perform the following checks, cycling the power after each.

Check all connectors are firmly pushed home. If any yellow is visible, then it is possible a good connection is not being made.

Check all cables for signs of visible damage. If a cable is thought to be damaged, then check it for continuity.

Finally, if the problem persists, then swap-out each cable in turn.

If the problem persists, contact your service center.

Inhibit Active

The system has detected that latched inhibit condition has occurred.

The origin of the inhibit depends on the wheelchair you are working with, but typical examples may include connection of a battery charger or the seat being in an elevated position.

Remove the likely source of the inhibit and cycle the power.

If the error persists and you are certain that the inhibit is not caused by an operating condition of the wheelchair, then it may be that an inhibit connector has become disconnected.

Inhibits may be connected to the Power Module (PM) or the Intelligent Seating/lighting Module (ISM). The PM contains Inhibits 2 and 3, while the ISM contains Inhibits 4 and 5. If you have access to the inhibit four digit error codes that start with "1E", then the last two digits represent the actual inhibit that is active. E.g. error code 1E03 relates to inhibit 3. If the last two digits are not 02, 03, 04 or 05, then the inhibit is probably being caused by a non-PGDT module.

Check the wiring to all the inhibits and cycle the power.

If the problem persists, contact your service center.

System Error

The system has detected an error, but cannot determine which module is reporting the error.

Check all cables and connections, and cycle the power.

If the problem persists and the system contains non-PGDT modules, then disconnect these modules and cycle the power.

If this clears the error, re-connect each non-PGDT module in turn until the error reappears. Once the defective module has been identified, contact your service center.

If the problem persists, then there is an error in one of the PGDT modules and you should contact your service center.

DIME Error

The system has detected some identification conflict between the modules that are connected.

If a new module has just been connected, disconnect it and cycle the power. If this clears the error, then there is a problem with the new module and you should contact your service center.

If the error persists, then disconnect each module in turn and cycle power. Once the defective module has been identified, contact your service center.

If it has not been possible to identify the module that is causing the error, then you should contact your service center.

Brake Fault Detect

Brake Fault Detect should only ever be turned off if there are no solenoid brakes fitted to the wheelchair.

Warning: PG Drives Technology accepts no liability for losses of any kind arising from failure to comply with this condition.

Brake Voltage

It is essential that the control system is programmed for the correct brake voltage. If it is not, damage may occur to the brakes or the drive performance of the wheelchair may be affected.

Soft Stop Rate

You must ensure that the emergency stopping distance is within the distance specified for the country in which the wheelchair will be used. TÜV Product Service (Germany) specify the distance to be as stated in EN12184.

Maximum Current Limit, Boost Drive Current, Boost Drive Time, Current Foldback Temperature, Minimum Current Limit

Each of these values has a maximum permissible value. These maximum values are shown in the table below and depend on the Power Module type you are working with. To establish the Power Module type, refer to the exact data sheet for the unit you are programming. These values should never be exceeded. Doing so will invalidate the Power Module's warranty and may affect the long term reliability of the system.

Parameter	PM60	PM80	PM120
Maximum Current Limit	60A	80A	100A
Boost Drive Current	60A	80A	120A
Boost Drive Time	-	-	10s
Current Foldback Temperature	60°C	60°C	60°C
Minimum Current Limit	50A	50A	50A

Compensation

The user of the software is responsible for ensuring that the controller is matched to the motor characteristics. Failure to do this may result in poor control characteristics, which in extreme instances can make a vehicle uncontrollable and potentially unsafe.

Never set Compensation to greater than 70% of the motor armature resistance. Failure to observe this warning could result in an unsafe set-up for the wheelchair user.

Warning: PG Drives Technology accepts no liability for losses of any kind arising from failure to comply with this condition.

Joystick Disconnected

The Joystick is disconnected. Please check connections and retry the operation.

Center Joystick

The most likely cause of this error is the joystick (or other type of Input Device) being away from center while the system is switched on. Ensure the joystick is in the center position when you switch the system on.

If the error persists, contact your service agent.

SID Disconnected

The Omni has detected that the Specialty Input Device (SID) has become disconnected. Check all cables and connectors between the Omni and the SID.

If the error persists, check that the setting of the parameter, 9-Way Detect, is appropriate for the SID that is being used. For example, if the SID has no detect-link, then this parameter should be set to Off.

If the error persists, contact your service agent.

User Switch Disconnected

The Omni has detected that the User Switch has become disconnected. Check all cables and connectors between the Omni and the User Switch.

If the error persists, contact your service agent.

If it is required to use the Omni without a User Switch being connected, then the parameter, Switch Detect, should be set to Off. If a User Switch is not used the responsibility for that decision lies with the healthcare professional.

Warning: Because a disconnected User Switch means there is no emergency stop function, PG Drives Technology recommend that Switch Detect is this option is always set to On. PG Drives Technology will accept no liability for any losses resulting from any other setting of this parameter.

Over Pressure

The Omni has detected that the pressure on the Sip and Puff transducer has exceeded its maximum level. Check the user is operating the Sip and Puff device correctly.

If the error persists, contact your service agent.

Cycle Power

Cycle the power to the system.

If the message persists, contact your service agent.

Language Error

One of the language files in the Power Module has become corrupted. Please contact your service agent.

No Trip Detected

No Trip Detected.

To detect a trip, the Dongle must already be connected when the trip occurs.

Switch off the system and repeat the sequence that led to the trip.

M1 Sensor Error

The control system has detected a problem connecting to the Encoder module via its M1 output. Whether the M1 output is for the left or the right motor depends on the configuration of the wheelchair and the programming of the Power Module.

Check all cables and connections between the Power Module M1 connector, the Encoder Module and the Encoder wiring, then cycle the power.

If the error persists, contact your service center.

M2 Sensor Error

The control system has detected a problem connecting to the Encoder module via its M2 output. Whether the M2 output is for the left or the right motor depends on the configuration of the wheelchair and the programming of the Power Module.

Check all cables and connections between the Power Module M2 connector, the Encoder Module and the Encoder wiring, then cycle the power.

If the error persists, contact your service center.

Encoder Module Disconnected

This occurs when the control system detects the Encoder Module has become disconnected. Check all cables and connections between the Power Module connector, the Encoder Module and the Encoder wiring, then cycle the power.

If the error persists, contact your service center.

Joystick Stationary Time and Joystick Stationary Range

The values of both these parameters have been selected by the wheelchair manufacturer in order to minimise the risk of motor damage. Do not adjust these values without consulting the wheelchair manufacturer.

Warning: PGDT accept no liability for losses of any kind resulting from the adjustment of these parameters.

Joystick Stationary Timeout

This occurs when the joystick is deemed to have been held stationary for an excessive period of time. The controller will stop drive to prevent possible damage the wheelchair's motors. Turning the control system Off and On again will clear this error message.

If the message persists contact your wheelchair manufacturer.

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